



MS STEYER
8-line hydrolicer



MECHANICAL TREATMENTS
Crew and Vessel

Highly qualified and skilled crew
with 4 years experience as a team.

Equipment mounted on board the
59-metre long and 15-metre wide MS Steyer.

Highest
standards
in delousing



TECHNICAL SPECIFICATION

VESSEL

Name: MS Steyer
Length: 59 meter
Width: 15 meter

HYGIENE

Surface washing: Central washing system
Foam / Disinfection / Flush

CIP wash: Hoses, pipes and pumps
Wash and Disinfection

Highest standards of hygiene with excellent
washing and disinfection procedures

OTHER

Catering: Hot lunch and dinner may be offered

CREW

10 crew members, including a chef
and a veterinarian/aqua medicin biologist

PROCESS SPECIFICATION

Hydrolicer model 6: 2 pcs per line
Number of lines: 8 pcs
Hydroflow: 3 pcs per line
Lifting height: 2.5-2.8 meters
Fish sizes: 0.5-15 kg
Fish transfer: Cage to cage or
back to source cage
(new feature, May 2022)

FILTER

Mechanical rotor filter: 3 pcs, 100 μ (micron)

CAPACITIES

Number of fish 1-3 kg: 60-70,000 fish/hour
Number of fish 3-12 kg: 55-60,000 fish/hour




AQUA PHARMA
A JOINT VENTURE OF AQUATIQ AND SOLVAY

For more information,
please visit our website:

www.aqua-pharma.com

Hovemoveien 1
2624 Lillehammer
NORWAY

Phone: +47 61 24 70 10
E-mail: post@aqua-pharma.com